

Faster Organochlorine Pesticide Sample Throughput

On New Rtx®-CLPesticides & Rtx®-CLPesticides2 Columns

By Jason Thomas, Environmental Innovations Chemist

- Dramatically improve sample throughput.
- Results in <7min. by conventional analysis, or <5min. using the Gerstel MACH system.
- Outstanding resolution on all columns.

As the environmental testing market continues to be very competitive, laboratory operating costs are a critical concern. Increasing sample throughput is one way to reduce costs, and shortening analytical run time is an effective way to do this. Here we offer methods for reducing run time for the organochlorine pesticides analyzed under US EPA Method 8081. The significant reduction in both analysis time and more significantly, cycle time, offered here is a major benefit for environmental laboratories.

Restek developed the Rtx®-CLPesticides and Rtx®-CLPesticides2 column pair specifically for chlorinated pesticides. These phases were designed to separate the isomers and the structurally similar pairs on the list of target analytes. Here we introduce new film thicknesses with optimized phase ratios for some of the columns in this line. Using these new stationary phase film thicknesses and the optimized run conditions shown, the 20 compounds in US EPA Method 8081 can be separated to baseline in less than 7 minutes (Figure 1). This allows rapid analysis without sacrificing column capacity, which translates, of course, into much improved sample throughput for your laboratory.

An Even Faster Alternative

In the attempt to obtain faster analytical run times, several different concepts have been introduced to improve the stock performance of standard GCs. One of the most recent and versatile ideas is the low thermal mass method by Gerstel using an apparatus called the MACH, (Modular Accelerated Column Heater) (Figure 2). This system operates by heating the capillary column outside of the GC oven in a small column module mounted on the oven door.

This apparatus provides several important advantages. First, due to the low thermal mass of the unit, very rapid heating and cooling times can be realized, which significantly shortens cycle times. Second, because of the way the column is wrapped, very uniform heating occurs, which eliminates the eddies and hot spots produced in a conventional GC oven. Finally, since the column modules are independently controlled, two different temperature programs can be run simultaneously, which allows each column to be optimized individually.

Restek applied this novel MACH technology to EPA Method 8081 using an Rtx®-CLPesticides and Rtx®-CLPesticides2 column pair. Almost 100% baseline resolution was obtained for all 22 pesticides and surrogates, on both columns, in under five minutes (Figure 3). This combination of ultra-fast analysis time and outstanding resolution is a result of the unique selectivity and high efficiency of the phases combined with the narrow peaks associated with ultra-rapid ramp rates.

Regardless of whether you choose to embrace the new fast-GC technology, or continue to adhere to more conventional GC, Restek Rtx®-CLPesticides and Rtx®-CLPesticides2 columns can provide exceptional performance and very rapid run times when analyzing chlorinated pesticides.






Figure 2 Two column modules in a Gerstel MACH column heating system.

The MACH system allows independent temperature programming of up to two columns, simultaneously.

thank you

Instrument provided courtesy of Gerstel USA.
www.gerstelus.com

Rtx[®]-CLPesticides Columns (fused silica)

ID	df (μm)	temp. limits	length	cat. #
0.18mm	0.18	-60 to 310/330°C	20-Meter	42102
0.53mm	0.50	-60 to 300/320°C	30-Meter	11140

Rtx[®]-CLPesticides2 Columns (fused silica)

ID	df (μm)	temp. limits	length	cat. #
0.18mm	0.14	-60 to 310/330°C	20-Meter	42302
0.53mm	0.42	-60 to 300/320°C	30-Meter	11340

Pesticide Surrogate Mix

decachlorobiphenyl 2,4,5,6-tetrachloro-*m*-xylene
200μg/mL each in acetone, 1mL/ampul
cat. # 32000

Organochlorine Pesticide Mix AB #2

(20 components)



aldrin	8μg/mL	dieldrin	16
α-BHC	8	endosulfan I	8
β-BHC	8	endosulfan II	16
δ-BHC	8	endosulfan sulfate	16
γ-BHC (lindane)	8	endrin	16
α-chlordane	8	endrin aldehyde	16
γ-chlordane	8	endrin ketone	16
4,4'-DDD	16	heptachlor	8
4,4'-DDE	16	heptachlor epoxide (B)	8
4,4'-DDT	16	methoxychlor	80

In hexane:toluene (1:1), 1mL/ampul
cat. # 32292

get connected

See page 20-21 for our list of connectors and connector kits.

Figure 3 Resolve organochlorine pesticides in less than 5 minutes using the Rtx[®]-CLPesticides columns and the MACH system.



Sample: Organochlorine Pesticide Mix AB #2 (8-80μg/mL each component in hexane:toluene 1:1, (cat.# 32292), Pesticide Surrogate Mix (200μg/mL each component in acetone, cat.# 32000)

Resprep™ Florisil® SPE Cartridges: Normal Phase

	3mL/500mg (50-pk.)	6mL/500mg (30-pk.)	6mL/1000mg (30-pk.)
Florisil®	24031	—	24034
(EPA SW 846 methods and CLP protocols)	24032*	26086**	26085**

*Teflon® frits **Glass tubes with Teflon® frits

CarboPrep™ SPE Cartridges

	Tube Volume, Bed Weight	qty.	cat.#
CarboPrep™ 90	3mL, 250mg	50-pk.	26091