



Assure LC/MS/MS System Performance for Drug Analyses

Using a System Suitability Test Mix

By Kristi Sellers, Clinical/Forensic Innovations Chemist and Houssain El Aribi, Ph.D., LC/MS Product Specialist, MDS ScieX

- Increase sample throughput and data quality with easy, reliable verification of LC/MS/MS performance.
- Extensively documented standard preparation assures accurate, consistent solutions.
- Method included in Cliquant® Drug Screen & Quant Software—automatically generates test reports.

Sample throughput is a critical issue in drug toxicology, and it can be adversely affected by inferior system performance. Poor system performance can produce unreliable data, increase downtime, and necessitate sample reanalysis, which ultimately decreases sample throughput. To ensure that your LC/MS/MS system is running properly, a system suitability mix should be analyzed on a regular basis before case samples are analyzed.


Restek and Applied Biosystems have developed a system suitability mix specifically for drug testing that contains compounds covering a wide range of molecular weights, polarities, and retention times (Table I). This standards mix is designed to verify system performance and identify system problems. Figure 1 shows a representative chromatogram (+MRM transitions) of this suitability mix analyzed on an Applied Biosystems API 3200™ LC/MS/MS system. This simple test evaluates the entire analytical system, including the autosampler, column, HPLC pumps, and mass spectrometer. The data is automatically compared to expected results by Applied Biosystem's Cliquant® Drug Screen & Quant Software to identify system problems.

Table I Mix components vary in chemical properties, providing a rigorous system performance test.

Mass Spectrometer Conditions:

Analyte	MW	RT (min)	Q1	Q3
Amiodarone	645	12.30	646.0	58.0
Amphetamine	135	4.21	136.1	91.1
Caffeine	194	1.72	195.1	122.9
Codeine	299	3.47	300.2	165.2
Diazepam	284	5.25	285.1	193.2
Doxepin	279	8.72	280.2	107.1
Haloperidol	375	9.08	376.1	123.0
Morphine	285	2.24	286.1	165.1

Figure 1 Increase sample throughput by verifying system readiness with a drug standard system suitability mix. (MRM transitions)



The Cliquid® Drug Screen & Quant Software automates this test and generates a verification report which highlights failures. Peak area, peak shape, retention time reproducibility, fragmentation, and library search function all are evaluated through the software by comparing the test mix data to expected results. For example, full scan linear ion trap MS/MS data for diazepam and caffeine are compared to the library to assess fragmentation. A mass spectral match of 80% or more must be achieved to pass this portion of the system suitability test. Otherwise, the failure will be highlighted on the automated report.

Use this system suitability mix for drug analyses to assure system performance and simplify troubleshooting.

Analyzing this system suitability mix for drug analysis on a regular basis assures system performance, improves data quality, increases sample throughput, and simplifies troubleshooting. Moreover, the Cliquid® Drug Screen & Quant Software for Routine Forensic Toxicology enables nonexpert LC/MS/MS users to employ this system suitability test with little effort.

Acknowledgement

Method and data supplied by Applied Biosystems.

References

H. El Arbi, T. Sasaki, A. Schreiber, K. Sellers, K. Herwehe. Development of an LC/MS/MS System Suitability Test for Forensic Toxicology Applications. Applied Biosystems/MDS Sciex, 2007.

Allure® PFP Propyl Columns (USP L43) Excellent Columns for LC/MS and ELSD

Physical Characteristics:

particle size: 5µm, spherical	endcap: fully endcapped
pore size: 60Å	pH range: 2.5 to 7.5
carbon load: 17%	temperature limit: 80°C

5µm Column, 2.1mm	cat. #	
30mm	9169532	
50mm	9169552	
5µm Column, 3.2mm	cat. #	
30mm	9169533	
50mm	9169553	
5µm Column, 2.1mm	cat. #	
30mm (with Trident Inlet Fitting)	9169532-700	
50mm (with Trident Inlet Fitting)	9169552-700	
5µm Column, 3.2mm	cat. #	
30mm (with Trident Inlet Fitting)	9169533-700	
50mm (with Trident Inlet Fitting)	9169553-700	
Allure® PFP Propyl Guard Cartridges	qty.	cat. #
10 x 2.1mm	3-pk.	916950212
10 x 4.0mm	3-pk.	916950210
20 x 2.1mm	2-pk.	916950222
20 x 4.0mm	2-pk.	916950220

ordering note

For other dimensions of these columns, visit our website at www.restek.com.

ABI/SCIEX Cliquid® Drug Screen Mix

Forensic Drug Screen Test Mixture

amiodarone	10µg/mL	diazepam	10
amphetamine	10	doxepine	10
caffeine	10	haloperidol	1
codeine	10	morphine	10

In P&T methanol, 1mL/ampul

cat. # 36340

Forensic Drug Screen Internal Standard

D5-diazepam D5-doxepine

10µg/mL each in P&T methanol, 10mL/ampul

cat. # 36341

Trident Direct Guard Cartridge System

Easy to Use, Low Dead Volume—The Ultimate Combination of Convenience and Column Protection



Trident Direct 20mm guard cartridge holder with filter

Protection against particulate matter and maximum protection against irreversibly adsorbed compounds.



Trident Direct 10mm guard cartridge holder with filter

Protection against particulate matter and moderate protection against irreversibly adsorbed compounds.

Description	qty.	cat.#
10mm guard cartridge holder with filter	ea.	25084
20mm guard cartridge holder with filter	ea.	25086
Connection tip for Waters-style end fittings	ea.	25088
PEEK tip standard fittings	ea.	25087

Get More!

Clinical/Forensics/Toxicology
Related Articles Online

“Fast Screening and Confirmation for Gamma-Hydroxybutyrate (GHB)”

www.restek.com/CFT

