



# Beyond C18—Increase Retention of Hydrophilic Compounds Using Biphenyl Columns

*Searching for a better way to retain hydrophilic aromatic drug compounds? Biphenyl phases, such as the Pinnacle® DB Biphenyl column, provide greater retention than alkyl phases. Use a Biphenyl column to separate difficult-to-retain polar aromatics from unretained matrix contaminants.*

By Amanda Rigdon, Pharmaceutical Innovations Chemist and Rick Lake, Pharmaceutical Market Development Manager


Many drug classes include compounds with aromatic ring structures, some of which also contain a sulfone or sulfoxide group. Both sulfur groups have dipole moments, adding a hydrophilic character to compounds containing these functional groups. The analysis of hydrophilic compounds on a traditional alkyl column (e.g., C18) can be problematic, since alkyl columns depend on hydrophobic (dispersive) interactions for retention. Since the sulfone and sulfoxide groups contain  $\pi$  bonds, the Biphenyl column's affinity toward compounds containing these bonds makes it a logical choice when increased retention of compounds containing these groups is desired.

To explore the selectivity of the biphenyl phase towards sulfur-containing aromatic compounds, phenyl sulfone, a simple probe, was analyzed on alkyl (C18), phenyl, phenyl hexyl, and Biphenyl columns to determine the relative retention of each phase, as measured by capacity factor ( $k'$ ). In order to ensure separation of analytes from unretained contaminants, a minimum  $k'$  value of 2 is recommended for most analyses, however in cases where there is little to no matrix interference, a  $k'$  of 1 may be acceptable. The data in Figure 1 show that phenyl sulfone is retained to a much greater degree on the Pinnacle® DB Biphenyl column, than on the other phases tested ( $k' = 2.08$ ). This is due to the unique retention mechanism of the biphenyl stationary phase, which can interact with both the hydrophobic aromatic ring and the hydrophilic sulfone group through  $\pi$ - $\pi$  interactions. Although the phenyl stationary phase also allows for the use of  $\pi$ - $\pi$  interactions, the biphenyl phase has a larger electron cloud and is significantly more retentive.

To further test the retention of the Biphenyl column, a second set of probes, consisting of compounds in the NSAID family, was analyzed. Tenoxicam, which contains a sulfone group, and sulfapyrazone, which contains a sulfoxide group, were analyzed along with a void marker (uracil). Although these compounds are more complex than the probe used in the first experiment, the same pattern of retention was observed (Figure 2). The Pinnacle® DB Biphenyl column exhibited the greatest retention for tenoxicam. With  $k'$  values of 0.33 on the C18 and 0.49 on the phenyl columns, tenoxicam shows almost no retention on these stationary phases. The phenyl hexyl phase performed slightly better with a  $k'$  value of 1.52 for tenoxicam. However, when tenoxicam was analyzed on the Biphenyl column under the same conditions, the  $k'$  value increased to 2.22, a value much more likely to provide adequate resolution from matrix components. Sulfapyrazone, a less polar compound, also followed the same pattern of retention (Table I).

The improved retention for hydrophilic aromatics shown here is due to the unique  $\pi$ - $\pi$  interaction retention mechanism of the Biphenyl phase. This mechanism is particularly useful for analysis of sulfone- and sulfoxide-containing drug compounds, which are not easily retained on alkyl or phenyl phases. The Biphenyl phase provides greater retention than alkyl and phenyl phases and is ideal for separating difficult-to-retain polar aromatics from unretained matrix contaminants.

**Figure 1** The Biphenyl phase is more retentive for phenyl sulfone than other alkyl and phenyl phases.





**Biphenyl columns are much more effective than alkyl, phenyl, or phenyl hexyl phases when increased retention of hydrophilic aromatics is desired.**

**Pinnacle® DB 1.9 $\mu$ m columns available!**  
[www.restek.com/uhplc](http://www.restek.com/uhplc)




**Figure 2** Only the Biphenyl phase retains both test probes to  $k' > 2$ , the level recommended to ensure separation from unrestrained matrix contaminants.




**A. Pinnacle® DB Biphenyl**


**More retention with  $< \frac{1}{2}$  the carbon load, compared to phenyl hexyl columns**




**B. phenyl**



**C. phenyl hexyl**



**D. C18**



**Table I** Biphenyl columns show improved retention of sulfone- and sulfoxide-containing aromatic drugs.

K' Value				
	Biphenyl	Phenyl hexyl	Phenyl	C18
Tenoxicam	2.23	1.39	0.637	0.235
Sulfipyrazone	4.18	3.90	1.88	1.89

**Product Listing**

**Pinnacle® DB Biphenyl Columns (USP L11)**

particle size: 1.9  $\mu$ m, 3  $\mu$ m or 5  $\mu$ m, spherical  
pore size: 140Å  
carbon load: 8%

endcap: yes  
pH range: 2.5 to 7.5  
temperature limit: 80°C

$\mu$ m Column, 1.0mm	cat. #	price
30mm	9409331	
50mm	9409351	
100mm	9409311	
150mm	9409361	
$\mu$ m Column, 2.1mm	cat. #	price
30mm	9409332	
50mm	9409352	
100mm	9409312	
150mm	9409362	
$\mu$ m Column, 3.2mm	cat. #	price
30mm	9409333	
50mm	9409353	
100mm	9409313	
150mm	9409363	
$\mu$ m Column, 4.6mm	cat. #	price
30mm	9409335	
50mm	9409355	
100mm	9409315	
150mm	9409365	
5 $\mu$ m Column, 1.0mm	cat. #	price
30mm	9409531	
50mm	9409551	
100mm	9409511	
150mm	9409561	
200mm	9409521	
250mm	9409571	
5 $\mu$ m Column, 2.1mm	cat. #	price
30mm	9409532	
50mm	9409552	
100mm	9409512	
150mm	9409562	
200mm	9409522	
250mm	9409572	
5 $\mu$ m Column, 3.2mm	cat. #	price
30mm	9409533	
50mm	9409553	
100mm	9409513	
150mm	9409563	
200mm	9409523	
250mm	9409573	
5 $\mu$ m Column, 4.6mm	cat. #	price
30mm	9409535	
50mm	9409555	
100mm	9409515	
150mm	9409565	
200mm	9409525	
250mm	9409575	

Also available in 1.9  $\mu$ m for UHPLC!  
Visit [www.restek.com](http://www.restek.com) for a complete listing.