

Rechargeable 'S' Traps

Molecular Sieve (1/4-inch: cat.# 20685, 1/8-inch: cat.# 20686)


Hydrocarbon (1/8-inch: cat.# 22016)

Reconditioning Instructions

To recondition the rechargeable trap:

1. Install the trap in the GC oven.
2. Set nitrogen flow between 50–70mL/min.
3. Condition at 300°C for 2–3 hours.

Note: Recondition every 6–12 months depending on carrier gas purity and usage.



**Call Technical Service at 800-356-1688 or 814-353-1300, ext. 4
(or your local Restek representative) if you have any questions
about this product or any other Restek product.**



#204-08-014

Rev. date: 11/05

Restek U.S. • 110 Benner Circle • Bellefonte, PA 16823

phone: 814-353-1300 • 800-356-1688 • fax: 814-353-1309 • www.restek.com

Restek France • phone: 33 (0)1 60 78 32 10 • fax: 33 (0)1 60 78 70 90 • e-mail: restekfr@club-internet.fr

Restek Ireland • phone: 44 2890 814576 • fax: 44 2890 814576 • e-mail: resteurope@aol.com

Thames Restek U.K. LTD • phone: 44 1494 563377 • fax: 44 1494 564990 • e-mail: sales@thamesrestek.co.uk

Restek GmbH • phone: +49 (0) 6172 2797 0 • fax: +49 (0) 6172 2797 77 • e-mail: info@restekgmbh.de