

Siltek®/Sulfinert® Deactivated Stainless Steel Tubing

Thank you for your purchase of our Siltek®/Sulfinert®-coated tubing. The Siltek®/Sulfinert coating is the most inert substrate available for the transfer of low levels of organosulfur compounds.

Instructions for Use

Maximum temperature: 360°C

Minimum suggested bending radius:

$1/16$ -inch OD = 1-inch radius

$1/8$ -inch OD = 2-inch radius

$1/4$ -inch OD = 4-inch radius

Handling: Handle the Siltek®/Sulfinert® tubing as you would any stainless steel tubing. Take precautions to avoid sharp bends or bends tighter than the minimum suggested bending radius. Any bend sharper than those listed may cause the tubing to stretch, potentially creating active sites as the coating layer density decreases. Additionally, we recommend passing a burst of nitrogen or clean compressed air through the tubing (at approximately 50 psi) to prevent damage from any particles released during the stress of bending. Once the tubing has been blown out and put into service, additional particle generation will not occur.

Cutting: Cut the tubing ends with a standard metal tubing cutter. Point the tubing end down when cutting or reaming to prevent metal filings from depositing inside the bore. Always blow out tubing sections after cutting and end-cleaning to ensure the removal of metal or coating particles generated during cutting.

Solvent compatibility: Siltek®/Sulfinert® tubing is compatible with a wide variety of organic solvents including methylene chloride, pentane, acetone, methanol, and water. Avoid hydrofluoric acid and bases or strong hydrochloric solutions. These acids and bases will damage the Siltek®/Sulfinert® layer.

**Call Technical Service at 800-356-1688 or 814-353-1300, ext. 4
(or your Restek representative) if you have any questions about this product or
any other Restek product.**



© 2009 Restek Corporation.

ISO 9001:2000 #302-02 [001]
Rev. date: 02/09

Restek U.S. • 110 Benner Circle • Bellefonte, PA 16823

phone: 814-353-1300 or 800-356-1688 • fax: 814-353-1309 • www.restek.com

Restek France • phone: 33 (0)1 60 78 32 10 • fax: 33 (0)1 60 78 70 90 • e-mail: restek@restekfrance.fr

Restek Ireland • phone: 44 2890 814576 • fax: 44 2890 814576 • e-mail: restekurope@aol.com

Thames Restek U.K. LTD • phone: 44 1494 563377 • fax: 44 1494 564990 • e-mail: sales@thamesrestek.co.uk

Restek GmbH • phone: +49 (0) 6172 2797 0 • fax: +49 (0) 6172 2797 77 • e-mail: info@restekgmbh.de