



Ultra II™ LC Columns

The Column Line Designed for Optimal Chromatography on **Any LC System**, Based on Highly Inert Restek Silica

Available Phases:

- C18
- Silica
- Biphenyl
- PFP Propyl
- Aromax
- Aqueous C18

Available Particle Sizes:

- 1.9µm for UHPLC
- 2.2µm for UFLC and RRLLC
- 3µm and 5µm for HPLC


**More phases
coming soon!**



Chromatography Products

www.restek.com

NEW! Ultra II™ LC Columns



Ultra II™ Aromax Columns (USP L11)

Ultra II Aromax is a unique reversed phase material that exhibits superior retention and selectivity for aromatic and/or unsaturated compounds, compared to conventional alkyl and phenyl phases. This column is a great alternative to our Biphenyl phase when increased retention is required. A very suitable choice for analysis of steroids, tetracyclines, drug metabolites, and other compounds that contain some degree of unsaturation.

Ultra II™ Biphenyl Columns (USP L11)

A unique reversed phase material that exhibits both increased retention and selectivity for aromatic and/or unsaturated compounds, compared to conventional alkyl and phenyl phases. This is a great alternative to a C18 column when alternative selectivity is desired. An excellent choice for the analysis of steroids, tetracyclines, drug metabolites, and other compounds that contain some degree of unsaturation.

Ultra II™ C18 Columns (USP L1)

A retentive, highly pure material that exhibits excellent peak shape for a wide range of compounds. This is a robust and very reproducible general-purpose reversed phase column.

Ultra II™ Aqueous C18 Columns (USP L1)

Highly retentive and selective for reversed phase separations of polar analytes. Highly base-deactivated. Compatible with highly aqueous (up to 100%) mobile phases.

Ultra II™ PFP Propyl Columns (USP L43)


A pentafluorophenyl phase with a propyl spacer. Highly retentive for basic analytes. An excellent phase for separating nucleosides, nucleotides, purines, pyrimidines, and halogenated compounds.

Ultra II™ Silica Columns (USP L3)


High surface area. Type B silica packing.

Ultra II™ is the first LC column line specifically designed for universal application—optimal chromatography on any system. This column line is built on a highly inert, high surface area silica that is completely Restek manufactured, providing excellent column-to-column reproducibility. Only Ultra II™ columns are available in a comprehensive range of particle sizes (1.9µm, 2.2µm, 3µm, and 5µm), creating truly scalable chromatography on any type of LC instrument, from conventional to ultra-high pressure systems. Ultra II™ columns are available in traditional phases (C18, Aqueous C18, Silica) and unique chemistries which provide alternate selectivity (Biphenyl, Aromax, PFP Propyl).

Parabens on Ultra II™ C18




Sulfa Drugs on Ultra II™ Biphenyl



Famotidine on Ultra II™ Aromax

Fast! 5µm

Peak List	Ret. Time (min.)
1. impurity A	3.779
2. famotidine	8.264
3. impurity D	9.180
4. impurity C	9.911
5. impurity B	14.018



Sample:
 Inj.: 10µL
 Conc.: 100µg/mL famotidine, 10µg/mL each impurity
 Sample diluent: methanol

Column: Ultra II™ Aromax
 Cat.#: 9607565
 Dimensions: 150mm x 4.6mm
 Particle size: 5µm
 Pore size: 100Å


Conditions:
 Instrument: Shimadzu Prominence UFLCXR
 Mobile phase: 20mM potassium phosphate (pH 2.5):methanol

Time (min.)	%B
0	5
15	35

Flow: 2.0mL/min.
 Temp.: 40°C
 Det.: UV @ 268nm

Faster! 3µm

Peak List	Ret. Time (min.)
1. impurity A	2.094
2. famotidine	4.460
3. impurity D	4.937
4. impurity C	5.334
5. impurity B	7.480



Sample:
 Inj.: 5µL
 Conc.: 100µg/mL famotidine, 10µg/mL each impurity
 Sample diluent: methanol

Column: Ultra II™ Aromax
 Cat.#: 9607313
 Dimensions: 100mm x 3.2mm
 Particle size: 3µm
 Pore size: 100Å


Conditions:
 Instrument: Shimadzu Prominence UFLCXR
 Mobile phase: 20mM potassium phosphate (pH 2.5):methanol

Time (min.)	%B
0	5
8	35

Flow: 1.2mL/min.
 Temp.: 40°C
 Det.: UV @ 268nm

Fastest! 2.2µm

Peak List	Ret. Time (min.)
1. impurity A	0.892
2. famotidine	2.082
3. impurity D	2.292
4. impurity C	2.517
5. impurity B	3.540



Sample:
 Inj.: 2µL
 Conc.: 100µg/mL famotidine, 10µg/mL each impurity
 Sample diluent: methanol

Column: Ultra II™ Aromax
 Cat.#: 9607853
 Dimensions: 50mm x 3.0mm
 Particle size: 2.2µm
 Pore size: 100Å


Conditions:
 Instrument: Shimadzu Prominence UFLCXR
 Mobile phase: 20mM potassium phosphate (pH 2.5):methanol

Time (min.)	%B
0	5
4	35

Flow: 1.2mL/min.
 Temp.: 40°C
 Det.: UV @ 268nm

Explosives on Ultra II™ Aromax

5µm



- | | |
|-------------------------------|---------------------------|
| 1. HMX | 10. PETN |
| 2. RDX | 11. 4-nitrotoluene |
| 3. nitroglycerin | 12. 1,3-dinitrobenzene |
| 4. nitrobenzene | 13. 2,6-dinitrotoluene |
| 5. 4-amino-2,6-dinitrotoluene | 14. tetryl |
| 6. 3,5-dinitroaniline | 15. 2,4-dinitrotoluene |
| 7. 2-nitrotoluene | 16. 1,3,5-trinitrobenzene |
| 8. 2-amino-4,6-dinitrotoluene | 17. 2,4,6-trinitrotoluene |
| 9. 3-nitrotoluene | |

LC_EV0484


Sample: Nitroaromatics and Nitramine Explosives by HPLC, EPA 8330B (cat.# 33204)
Inj.: 10µL
Conc.: 10µg/mL each component
Sample diluent: methanol

Conditions:
Instrument: Shimadzu Prominence UFLCXR
Mobile phase: water:methanol, 35:65 (v/v)
Flow: 1.2mL/min.
Temp.: 30°C
Det.: UV @ 254nm and 210nm

Column: Ultra II™ Aromax
Cat.#: 9607575
Dimensions: 250mm x 4.6mm
Particle size: 5µm
Pore size: 100Å

Fat Soluble Vitamins on Ultra II™ C18

2.2µm



- | Peak List | Ret. Time (min.) |
|----------------------|------------------|
| 1. vitamin K2 | 0.905 |
| 2. vitamin D2 | 1.224 |
| 3. vitamin D3 | 1.305 |
| 4. vitamin E acetate | 1.828 |
| 5. vitamin K1 | 2.169 |

LC_PH0492

Sample:
Inj.: 1µL
Conc.: 100µg/mL each component
Sample diluent: acetone

Conditions:
Instrument: Shimadzu Prominence UFLCXR
Mobile phase: acetonitrile:methanol (85:15)
Flow: 1.5mL/min.
Temp.: ambient
Det.: UV @ 230nm


Column: Ultra II™ C18
Cat.#: 9604853
Dimensions: 50mm x 3.0mm
Particle size: 2.2µm
Pore size: 100Å

NSAIDs on Ultra II™ Biphenyl

2.2μm

Peak List

Peak	Ret. Time (min.)
1. ketoprofen	0.894
2. piroxicam	1.068
3. sulindac	1.197
4. diclofenac	1.465



LC_PH0487

Sample:

Inj.: 5μL
 Conc.: 15μg/mL each component
 Sample diluent: methanol

Column:

Cat.#: Ultra II™ Biphenyl
 9609853
 Dimensions: 50mm x 3.0mm
 Particle size: 2.2μm
 Pore size: 100Å

Conditions:

Instrument: Shimadzu Prominence UFLCXR
 Mobile phase: 20mM potassium phosphate,
 pH 2.5 with H₃PO₄:methanol (30:70)


Flow: 1.2mL/min.
 Temp.: 40°C
 Det.: UV @ 254nm

Organic Acids on Ultra II™ Aqueous C18

5μm

Peak List

Peak	Ret. Time (min.)
1. tartaric acid	2.061
2. quinic acid	2.210
3. malic acid	2.596
4. citric acid	4.177
5. fumaric acid	5.853



LC_PH0498

Sample:

Inj.: 5μL
 Conc.: 10μg/mL fumaric acid,
 2,000μg/mL each other acids
 Sample diluent: water

Column:

Cat.#: Ultra II™ Aqueous C18
 9608565
 Dimensions: 150mm x 4.6mm
 Particle size: 5μm
 Pore size: 100Å

Conditions:

Instrument: Shimadzu Prominence UFLCXR
 Mobile phase: 100% 20mM potassium phosphate (pH 2.5)
 Flow: 1.0mL/min.
 Temp.: 30°C
 Det.: UV @ 226nm

Alcohol Metabolites on Ultra II™ Biphenyl

5µm

Sample:
 Inj.: 5µL
 Conc.: 5µg/mL each component
 Sample diluent: urine diluted 1:10 with mobile phase


Column: Ultra II™ Biphenyl
 Cat.#: 9609552
 Dimensions: 50mm x 2.1mm
 Particle size: 5µm
 Pore size: 100Å

Conditions:
 Instrument: Shimadzu Prominence UFLCXR
 Mobile phase: A: 5mM dihexylammonium acetate in water
 B: 5mM dihexylammonium acetate in methanol

Time (min.)	%B
0.0	40
2.5	40
2.6	100
5.4	100
5.5	40
7.5	40

Flow: 0.3mL/min.
 Temp.: 40°C
 Det.: Shimadzu LCMS-2010EV
 Ion source: EST-
 Mode: SIM
 Source temp.: 200°C
 Source voltage: -4,000V

Peak List	m/z
1. matrix	125
2. ethyl glucuronide	221
3. ethyl sulfate	125



Cannabinoids on Ultra II™ Biphenyl

2.2µm

Sample:
 Inj.: 20µL
 Conc.: 300ng/mL urine sample with 50ng/mL internal standard

Column: Ultra II™ Biphenyl
 Cat.#: 9609853
 Dimensions: 50mm x 3.0mm
 Particle size: 2.2µm
 Pore size: 100Å

Conditions:
 Instrument: Shimadzu Prominence UFLCXR
 Mobile phase: A: 0.1% formic acid in water
 B: 0.1% formic acid in methanol

Time (min.)	%B
0.0	50
0.5	50
3.5	100
4.5	100
4.6	50
6.0	50


Flow: 0.7mL/min.
 Temp.: 35°C
 Det.: Applied Biosystems 3200 QTRAP™ LC/MS/MS system
 Ion source: TurboIonSpray®, ESI+
 Mode: MRM
 Source temp.: 550°C
 Source voltage: 5,000V
 Dwell: 50ms

Peak List:

- cannabidiol
- THC-COOH
- cannabinol
- Δ⁹-THC-d₃
- Δ⁹-THC

Mass Spectrometer Experiments:

Analyte	Q1	Q3
cannabidiol 1	315.2	123.1
cannabidiol 2	315.2	193.2
THC-COOH 1	345.2	299.3
THC-COOH 2	345.2	193.2
cannabinol 1	311.2	223.2
cannabinol 2	311.2	178.2
Δ ⁹ -THC-d ₃ 1	318.2	196.2
Δ ⁹ -THC 1	315.2	193.2
Δ ⁹ -THC 2	315.2	123.1



Ultra II™ Silica Columns (USP L3)

Physical Characteristics:

particle size: 1.9µm, 2.2µm, 3µm or 5µm, spherical
pore size: 100Å
carbon load: 0%

endcap: no
pH range: 2.5 to 7.5
temperature limit: 80°C

1.9µm Column, 2.1mm	cat. #
30mm	9600232
50mm	9600252
100mm	9600212
2.2µm Column, 3.0mm	cat. #
30mm	9600833
50mm	9600853
100mm	9600813
3µm Column, 1.0mm	cat. #
30mm	9600331
50mm	9600351
100mm	9600311
150mm	9600361
3µm Column, 2.1mm	cat. #
30mm	9600332
50mm	9600352
100mm	9600312
150mm	9600362
3µm Column, 3.2mm	cat. #
30mm	9600333
50mm	9600353
100mm	9600313
150mm	9600363
3µm Column, 4.6mm	cat. #
30mm	9600335
50mm	9600355
100mm	9600315
150mm	9600365
5µm Column, 1.0mm	cat. #
30mm	9600531
50mm	9600551
100mm	9600511
150mm	9600561
200mm	9600521
250mm	9600571
5µm Column, 2.1mm	cat. #
30mm	9600532
50mm	9600552
100mm	9600512
150mm	9600562
200mm	9600522
250mm	9600572
5µm Column, 3.2mm	cat. #
30mm	9600533
50mm	9600553
100mm	9600513
150mm	9600563
200mm	9600523
250mm	9600573
5µm Column, 4.6mm	cat. #
30mm	9600535
50mm	9600555
100mm	9600515
150mm	9600565
200mm	9600525
250mm	9600575

Ultra II™ Biphenyl Columns (USP L11)

Physical Characteristics:

particle size: 2.2µm, 3µm or 5µm, spherical
pore size: 100Å
carbon load: 15%

endcap: fully endcapped
pH range: 2.5 to 7.5
temperature limit: 80°C

2.2µm Column, 3.0mm	cat. #
30mm	9609833
50mm	9609853
100mm	9609813
3µm Column, 1.0mm	cat. #
30mm	9609331
50mm	9609351
100mm	9609311
150mm	9609361
3µm Column, 2.1mm	cat. #
30mm	9609332
50mm	9609352
100mm	9609312
150mm	9609362
3µm Column, 3.2mm	cat. #
30mm	9609333
50mm	9609353
100mm	9609313
150mm	9609363
3µm Column, 4.6mm	cat. #
30mm	9609335
50mm	9609355
100mm	9609315
150mm	9609365
5µm Column, 1.0mm	cat. #
30mm	9609531
50mm	9609551
100mm	9609511
150mm	9609561
200mm	9609521
250mm	9609571
5µm Column, 2.1mm	cat. #
30mm	9609532
50mm	9609552
100mm	9609512
150mm	9609562
200mm	9609522
250mm	9609572
5µm Column, 3.2mm	cat. #
30mm	9609533
50mm	9609553
100mm	9609513
150mm	9609563
200mm	9609523
250mm	9609573
5µm Column, 4.6mm	cat. #
30mm	9609535
50mm	9609555
100mm	9609515
150mm	9609565
200mm	9609525
250mm	9609575

1.9µm coming soon!

Ultra II™ Aromax Columns (USP L11)

Physical Characteristics:

particle size: 3µm or 5µm, spherical
pore size: 100Å
carbon load: 17%

endcap: fully endcapped
pH range: 2.5 to 7.5
temperature limit: 80°C

3µm Column, 1.0mm	cat. #
30mm	9607331
50mm	9607351
100mm	9607311
150mm	9607361
3µm Column, 2.1mm	cat. #
30mm	9607332
50mm	9607352
100mm	9607312
150mm	9607362
3µm Column, 3.2mm	cat. #
30mm	9607333
50mm	9607353
100mm	9607313
150mm	9607363
3µm Column, 4.6mm	cat. #
30mm	9607335
50mm	9607355
100mm	9607315
150mm	9607365
5µm Column, 1.0mm	cat. #
30mm	9607531
50mm	9607551
100mm	9607511
150mm	9607561
200mm	9607521
250mm	9607571
5µm Column, 2.1mm	cat. #
30mm	9607532
50mm	9607552
100mm	9607512
150mm	9607562
200mm	9607522
250mm	9607572
5µm Column, 3.2mm	cat. #
30mm	9607533
50mm	9607553
100mm	9607513
150mm	9607563
200mm	9607523
250mm	9607573
5µm Column, 4.6mm	cat. #
30mm	9607535
50mm	9607555
100mm	9607515
150mm	9607565
200mm	9607525
250mm	9607575

1.9µm & 2.2µm coming soon!

ordering note

For guard cartridges for these columns, visit our website at www.restek.com.

Visit us at www.restek.com/ultra2 for our most complete listing of Ultra II™ columns.

Ultra II™ C18 Columns (USP L1)

Physical Characteristics:

particle size: 1.9µm, 2.2µm, 3µm or 5µm, spherical
pore size: 100Å
carbon load: 19%
endcap: fully endcapped
pH range: 2.5 to 7.5
temperature limit: 80°C

1.9µm Column, 2.1mm	cat. #
30mm	9604232
50mm	9604252
100mm	9604212
2.2µm Column, 3.0mm	cat. #
30mm	9604833
50mm	9604853
100mm	9604813
3µm Column, 1.0mm	cat. #
30mm	9604331
50mm	9604351
100mm	9604311
150mm	9604361
3µm Column, 2.1mm	cat. #
30mm	9604332
50mm	9604352
100mm	9604312
150mm	9604362
3µm Column, 3.2mm	cat. #
30mm	9604333
50mm	9604353
100mm	9604313
150mm	9604363
3µm Column, 4.6mm	cat. #
30mm	9604335
50mm	9604355
100mm	9604315
150mm	9604365
5µm Column, 1.0mm	cat. #
30mm	9604531
50mm	9604551
100mm	9604511
150mm	9604561
200mm	9604521
250mm	9604571
5µm Column, 2.1mm	cat. #
30mm	9604532
50mm	9604552
100mm	9604512
150mm	9604562
200mm	9604522
250mm	9604572
5µm Column, 3.2mm	cat. #
30mm	9604533
50mm	9604553
100mm	9604513
150mm	9604563
200mm	9604523
250mm	9604573
5µm Column, 4.6mm	cat. #
30mm	9604535
50mm	9604555
100mm	9604515
150mm	9604565
200mm	9604525
250mm	9604575

Ultra II™ Aqueous C18 Columns (USP L1)

Physical Characteristics:

particle size: 2.2µm, 3µm or 5µm, spherical
pore size: 100Å
carbon load: 15%
endcap: no
pH range: 2.5 to 7.5
temperature limit: 80°C

2.2µm Column, 3.0mm	cat. #
30mm	9608833
50mm	9608853
100mm	9608813
3µm Column, 1.0mm	cat. #
30mm	9608331
50mm	9608351
100mm	9608311
150mm	9608361
3µm Column, 2.1mm	cat. #
30mm	9608332
50mm	9608352
100mm	9608312
150mm	9608362
3µm Column, 3.2mm	cat. #
30mm	9608333
50mm	9608353
100mm	9608313
150mm	9608363
3µm Column, 4.6mm	cat. #
30mm	9608335
50mm	9608355
100mm	9608315
150mm	9608365
5µm Column, 1.0mm	cat. #
30mm	9608531
50mm	9608551
100mm	9608511
150mm	9608561
200mm	9608521
250mm	9608571
5µm Column, 2.1mm	cat. #
30mm	9608532
50mm	9608552
100mm	9608512
150mm	9608562
200mm	9608522
250mm	9608572
5µm Column, 3.2mm	cat. #
30mm	9608533
50mm	9608553
100mm	9608513
150mm	9608563
200mm	9608523
250mm	9608573
5µm Column, 4.6mm	cat. #
30mm	9608535
50mm	9608555
100mm	9608515
150mm	9608565
200mm	9608525
250mm	9608575

1.9µm coming soon!

PATENTS & TRADEMARKS

Restek patents and trademarks are the property of Restek Corporation. Other trademarks appearing in Restek literature or on its website are the property of their respective owners.

Ultra II™ PFP Propyl Columns (USP L43)

Physical Characteristics:

particle size: 2.2µm, 3µm or 5µm, spherical
pore size: 100Å
carbon load: 11%
endcap: fully endcapped
pH range: 2.5 to 7.5
temperature limit: 80°C

2.2µm Column, 3.0mm	cat. #
30mm	9606833
50mm	9606853
100mm	9606813
3µm Column, 1.0mm	cat. #
30mm	9606331
50mm	9606351
100mm	9606311
150mm	9606361
3µm Column, 2.1mm	cat. #
30mm	9606332
50mm	9606352
100mm	9606312
150mm	9606362
3µm Column, 3.2mm	cat. #
30mm	9606333
50mm	9606353
100mm	9606313
150mm	9606363
3µm Column, 4.6mm	cat. #
30mm	9606335
50mm	9606355
100mm	9606315
150mm	9606365
5µm Column, 1.0mm	cat. #
30mm	9606531
50mm	9606551
100mm	9606511
150mm	9606561
200mm	9606521
250mm	9606571
5µm Column, 2.1mm	cat. #
30mm	9606532
50mm	9606552
100mm	9606512
150mm	9606562
200mm	9606522
250mm	9606572
5µm Column, 3.2mm	cat. #
30mm	9606533
50mm	9606553
100mm	9606513
150mm	9606563
200mm	9606523
250mm	9606573
5µm Column, 4.6mm	cat. #
30mm	9606535
50mm	9606555
100mm	9606515
150mm	9606565
200mm	9606525
250mm	9606575

1.9µm coming soon!



Lit. Cat.# GNTS1177

© 2009 Restek Corporation.

Restek U.S. • 110 Benner Circle • Bellefonte, PA 16823 • 814-353-1300 • 800-356-1688 • fax: 814-353-1309 • www.restek.com

Restek France • phone: +33 (0)1 60 78 32 10 • fax: +33 (0)1 60 78 70 90 • e-mail: restek@restekfrance.fr

Restek GmbH • phone: +49 (0)6172 2797 0 • fax: +49 (0)6172 2797 77 • e-mail: info@restekgmbh.de

Restek Ireland • phone: +44 (0)2890 814576 • fax: +44 (0)2890 814576 • e-mail: restekeuropa@aol.com

Restek Japan • phone: +81 (3)6459 0025 • fax: +81 (3)6459 0025 • e-mail: ryosei.kanaguchi@restek.com

Thames Restek U.K. LTD • phone: +44 (0)1494 563377 • fax: +44 (0)1494 564990 • e-mail: sales@thamesrestek.co.uk

