

Overcoming the Challenges of Hand Sanitizer Analysis

Chris English¹, Tyler West², Bobby Polak², Michael Ratkovich², Brian Sloat², and Michael Sandova².
¹ Restek Corporation, ² Santé Laboratories

Introduction

The U.S. Food and Drug Administration (FDA) has placed methanol containing hand sanitizers on an import alert. These products do not list this ingredient and some are incorrectly labeled "FDA approved." Nearly 100 hand sanitizers have been found to contain methanol. Methanol, known as wood alcohol, is commonly found in solvents, and windshield washer fluid. Ingestion of 10 mL can cause permanent blindness and 15 mL is considered a lethal dose. Transdermal methanol poisoning has been well documented and may result in optic nerve necrosis with permanent eye damage. Denatured alcohol is ethanol mixed with other alcohols and can contain 50% methanol, along with dozens of other possible contaminants.

The separation and analysis of these active and inactive compounds is challenging and finding the correct column and conditions to implement and validate a robust test method is difficult for any method development team. We started with computer software capable of modeling these compounds on a variety of stationary phases which allows us to adjust; flow, temperature, film thickness, column diameter and column length. Using the column recommendation, Santé Laboratories was able to transfer, qualify, optimize and validate a GC-MS published FDA method. The method parameters, chromatograms, and challenges encountered in developing this method, will be presented in addition to an introduction to novel modeling software.

Using Computer Modeling Software to Simulate Gas Chromatographic Separations

A free to use, web-based computer program has been employed to accelerate Gas Chromatography (GC) method development by direct simulation of the chromatographic process. The program makes use of a highly accurate time-summation modeling approach, coupled with large libraries of data to provide absolute retention time predictions.

The user is able to directly control every parameter of the chromatography being modeled including carrier gas, stationary phase, column dimensions, and temperature programs. The temperature program can be iteratively optimized to obtain the desired separation in the minimum run time. Pro EZGC is preloaded for modeling of more than 20 different stationary phases and thousands of compounds including; Solvents, Pesticides, PCB's, PAH's, semivolatiles, and volatile organics.

Pro EZGC® Chromatogram Modeler

Compounds Conditions My EZGC®

Search by Name or CAS # Search by Phase >>

acetaldehyde
 methanol ← #1 Enter compounds here
 ethanol
 acetone
 2-propanol
 1-propanol
 ethyl acetate
 2-butanol

Clear Solve ← #2 Solve

Results were found on 3 stationary phases:

Rxi-624Sil MS (14 out of 14 resolved)
 Rxi-624Sil MS (14 out of 14 resolved) ← Columns available
 Stabilwax (14 out of 14 resolved)
 Rtx-VMS (14 out of 14 resolved)

Figure 1: Using the FDA's analytical procedure to assess finished hand sanitizers; Level 1 & Level II compounds were entered into Pro EZGC for modeling. Notice three columns were able to baseline resolve these compounds: Rxi-624Sil MS, Rtx-VMS and the Stabilwax.




Figure 2: Pro EZGC Model solved for FDA's Level I & Level II compounds where chromatographic conditions were entered, including; carrier gas flow rate, vacuum, and oven temperature gradient.




Figure 3: Chromatogram generated by Santé Laboratories using the FDA analytical procedure to assess finished hand sanitizers; Level 1 & Level II compounds closely match the Pro EZGC model (figure 2).

Conditions and Instrumentation

Instrument: Thermo Scientific Trace 1310 GC, ISQ 7000 MS
 TriPlus RSH Autosampler
 Column: Rxi-624Sil MS, 30.00 m, 0.25 mm ID, 1.40 µm
 Carrier Gas: Helium, Constant Flow @ 1.00 mL/min, 36.12 cm/sec
 Temperature Gradient: 40 °C (hold 5 min) to 240 °C @ 30 °C/min (hold 4 min)
 Injection: 1µL, Split 50:1, 250°C
 Liner: Topaz Precision Split Liner with Wool

Direct Injection GC-MS Method for the Detection of Listed Impurities in Hand Sanitizers

The FDA's analytical procedure can be used to evaluate products with either ethanol or isopropyl alcohol labeled as the active ingredient to assess percent alcohol as well as target impurities. This approach utilized MS, a departure from the USP <467> which allows detection by FID. In addition spike and recovery studies are conducted to evaluate matrix effects and standards are prepared using a stock standard working solution in acetonitrile.

Calibration Stocks	Mixed Stock		Working Standard (Nominal)	Spiked Blank (Nominal)	Working Standard (Measured)						AVG	STDEV	RSD	
	mg	mL	ug/mL	ug/mL	WS1	WS2	WS3	WS4	WS5	WS6				
1,1-Diethoxyethane (Acetal)	10.91	100	109.1	10.9	1.1	11.3	10.9	10.5	10.8	10.8	9.4	10.6	0.6	5.4
1-Butanol	152.63	100	1526.3	152.6	15.3	158.3	151.3	148.1	149.9	151.0	132.5	148.5	7.8	5.3
1-Pentanol (n-Amyl Alcohol)	147.02	100	1470.2	147.0	14.7	155.0	141.7	138.0	146.6	146.9	129.9	143.0	7.9	5.5
1-Propanol	155.15	100	1551.5	155.2	15.5	160.3	155.0	151.0	154.7	154.3	137.5	152.1	7.1	4.6
2-Butanol	144.80	100	1448.0	144.8	14.5	151.3	144.7	140.8	143.0	142.4	122.7	140.8	8.7	6.2
2-Methyl-1-Propanol (Isobutanol)	146.10	100	1461.0	146.1	14.6	152.6	144.0	141.9	144.5	143.2	125.4	141.9	8.2	5.8
2-Propanol (Isopropanol)	380.54	100	3805.4	380.5	38.1	398.8	376.1	366.7	373.9	368.0	321.7	367.5	23.1	6.3
3-Methyl-1-Butanol	158.69	100	1586.9	158.7	15.9	164.2	157.4	154.5	157.4	154.7	137.1	154.2	8.3	5.4
Acetaldehyde	10.49	100	104.9	10.5	1.0	11.0	10.2	10.1	10.4	10.1	8.4	10.1	0.8	8.0
Acetone	154.89	100	1548.9	154.9	15.5	163.8	153.2	149.1	151.1	149.0	128.0	149.0	10.7	7.2
Benzene, ~1000 ug/mL in ACN	392.58	100	3411.4	0.3	0.0	0.4	0.3	0.3	0.3	0.3	0.3	0.3	0.0	5.5
Ethanol	388.27	100	3882.7	388.3	38.8	402.7	383.2	379.7	387.7	385.1	336.8	379.2	20.3	5.4
Ethyl Acetate	173.77	100	1737.7	173.8	17.4	182.4	172.5	168.9	170.6	169.7	148.1	168.7	10.3	6.1
Methanol	152.51	100	1525.1	152.5	15.3	159.3	151.7	147.8	152.4	150.4	130.9	148.8	8.7	5.9

Figure 4: Quantitation of the impurities is based on peak area of impurity in sample chromatogram compared to an external calibration standard containing the reference standards. These working standards were analyzed for standard deviation and relative standard deviations (%RSD) where the replicates may not exceed 10% RSD.




Figure 5: Client sample of finished product exceeding acceptable levels of acetaldehyde (86.6 ug/g) and below limit concentrations of methanol (less than 316.4 ug/ml). Other impurities match retention times and spectra but fall well below FDA limit guidelines.

Conclusions

The compliance testing team at Santé Laboratories; a contract development and manufacturing organization (CDMO) headquartered in Austin, Texas, performs raw materials and finished product lot release testing for a variety of industries. The testing includes active ingredients and excipients frequently used to manufacture effective and safe hand sanitizers. Restek's Pro EZGC software enables laboratories to quickly determine GC conditions for targeted methods and specific resolution requirements. This software was successfully used to predict retention times using FDA's GC conditions and Santé Laboratories was able to develop and validate this analytical procedure.