

# Genuine Restek Replacement Parts for Beckman and Hitachi HPLC Systems

Keep Your System Running Smoothly

by Rebecca Wittrig, Ph.D., HPLC Product Marketing Manager

- ✓ Designed to meet or exceed original equipment performance.
- ✓ Simplify ordering—a single source for columns, tools, parts, reference mixes.
- ✓ Renowned Plus 1™ service.



## Restek Replacement Parts for Beckman HPLC Systems

Description	Model #	Similar to Beckman part #	qty.	cat.#	price
Inlet Check Valve Assembly	110A&B, 112, 114M, 116, 118, 125, 126, 127, 128	240720	ea.	25439	
Outlet Check Valve Assembly	110A&B, 112, 114M, 116, 118, 125, 126, 127, 128	240721	ea.	25440	
Inlet Check Valve Cartridge for Check Valve Assembly	110A&B, 112, 114M, 116, 118, 125, 126, 127, 128	240620	ea.	25441	
Outlet Check Valve Cartridge for Check Valve Assembly	110A&B, 112, 114M, 116, 118, 125, 126, 127, 128	240621	ea.	25442	
Graphite Guide (Bushing)	Pumps	243714	ea.	25443	
Guide Sleeve for Graphite Plunger Guide	110 Series	243713	ea.	25444	
Outlet Check Filter Frit Assembly	Pumps	240619	ea.	25445	
Piston Guide Assembly	Pumps	243045	ea.	25446	
Plunger	110 Series and 112 Pumps	243053	ea.	25447	
Plunger Seal	110 Series	887138	ea.	25448	
Plunger Seal	112 Pumps	236797	ea.	25449	
Sapphire Plunger	114M, 116, 118, 125, 126, 127, 128 Pumps	240714	ea.	25450	
Plunger Seal	114M, 116, 118, 125, 126, 127, 128 Pumps	241037	ea.	25451	
Pump Seal, Gold	114, 116, 125, 126, 127, 128 Pumps	241037	ea.	25452	
Plunger Wash Seal	Pumps	238627	ea.	25453	
Deuterium Lamp	DU60, 62, 64, 65	596791	ea.	25454	

## Restek Replacement Parts for Hitachi HPLC Systems

Description	Model #	Similar to Hitachi part #	qty.	cat.#	price
SS Check Valve Cartridges (1/16" Ruby Ball)	655/6000/6200	ANO-085	2-pk.	25455	
Inlet Check Valve Assembly	655/6000/6200	ANO-0833	ea.	25456	
Outlet Check Valve Assembly	655/6000/6200	ANO-0834	ea.	25457	
Plunger Assembly	655, 6000, 6200, 7100	810-1033	ea.	25458	
Plunger Assembly	L655A, 6000, 6200, 7100	655-1080	ea.	25459	
Plunger Assembly, Gold	655, 6000, 6200, 7100	655-1080	ea.	25460	
Inlet Check Valve Assembly	L-7100	ANO-0836	ea.	25461	
Outlet Check Valve Assembly	L-7100	ANO-0837	ea.	25462	
Rotor Assembly for Dilutor Valve	AS-7200 Autosampler	810-3085	ea.	25463	
Rotor Seal Kit	AS-7200, AS-7250 Injection Valve	ANO-0818	ea.	25464	
Deuterium Lamp, Prealigned	L4000, L4200, L4250, L7400	885-3570	ea.	25465	



### Sonic Debubbler

- ✓ Fast and neat—no breaking connections.
- ✓ Easy to use.
- ✓ Less solvent waste; less clean-up.

**COOL TOOL!**

Just touch the Sonic Debubbler to the inlet line or check valve — sonic vibrations will quickly dislodge or redissolve trapped air bubbles. Reduces downtime or conversion time from one mobile phase to another.

Description	qty.	cat.#	price
Sonic Debubbler	ea.	20444	
Sonic Debubbler (220V)	ea.	25098	



### Genuine Restek Replacement Parts for HPLC Systems

(lit. cat.# 59012)  
A listing of Restek parts for Agilent, Beckman, Hitachi, PerkinElmer, Shimadzu, and Waters instruments.



### HPLC Column Selection Guide

(lit. cat.# 59454A)  
A useful chart to keep with your workbooks, or post on a wall. Quickly scan important characteristics of Restek HPLC columns. Includes cross-references to similar phases.



### HPLC Tech Tips Wall Chart

(lit. cat.# 59894A)  
Almost everything you need to remember about HPLC, condensed into 3 feet by 2 feet: mobile phase basics, buffers (types, pK<sub>a</sub>s, pH ranges, formula masses, more), miscibility and solubility chart (invaluable!), system setup and optimization,

detector tips, pressure conversion factors, most-used chromatographic equations, column storage essentials. Post near your instrument to save time; perhaps save a column.

**Free!**

Call our literature hotline: 800-356-1688 or 814-353-1300, ext. 5, or contact your Restek representative.