

new! **HPLC Autosampler Syringes**
Now Available from Restek—Convenient, One-Stop Shopping

by Donna Lidgett, GC Accessories Product Marketing Manager

- Made by leading syringe manufacturers.
- Designed and tested for reliable performance.
- Fast delivery.

These syringes join an extensive selection of GC and HPLC syringes available from Restek. To see other syringes, refer to our catalog, or visit our website.



free literature

Genuine Restek
Replacement Parts for
HPLC Systems

Defective parts can cause poor chromatography and shorten column lifetimes. Use Restek check valves, pump piston seals, detector lamps, and other parts to keep your Agilent, Beckman, Hitachi, PerkinElmer, Shimadzu, or Waters system running smoothly and chromatographically sharp. Genuine Restek Replacement Parts equal or exceed the performance of original equipment components.
Request lit. cat.# 59012 (4 pages).

SGE Syringes for Agilent 1090 & 1100 LC Autosamplers

- Termination: 1/4-32 UNF thread



Volume	SGE			Restek	
	Model	cat.#	qty.	cat.#	
25µL	25D-HP1090-GT	003670	ea.	22290	
250µL	250D-HP1090-GT	006670	ea.	22291	

SGE Syringe for Hitachi LC Autosamplers

- Termination: M10x1 thread



Volume	SGE			Restek	
	Model	cat.#	qty.	cat.#	
500µL	500C-HITACHI	007660	ea.	22292	

SGE Syringes for PerkinElmer LC Autosamplers

- Termination: 1/4-28 UNF thread



Volume	SGE			Restek	
	Model	cat.#	qty.	cat.#	
50µL	50D-CX-GT	004995	ea.	22295	
100µL	100D-CX-GT	005990	ea.	22296	
250µL	250D-CX-GT	006995	ea.	22297	
500µL	500D-CX-GT	007995	ea.	22298	
1mL	1MD-C-GT	008185	ea.	22299	

SGE Syringes for Waters™ WISP® LC Autosamplers

- Termination: 1/4-28 UNF thread



Volume	SGE			Restek	
	Model	cat.#	qty.	cat.#	
25µL	25D-WISP	003990	ea.	22293	
250µL	250D-WISP	006690	ea.	22294	

Hamilton Syringes for Waters™ WISP® LC Autosamplers

- Teflon® PTFE-tipped plunger

Volume	Hamilton			Restek	
	Model	cat.#	qty.	cat.#	
25µL	1702	80020	ea.	24528	
250µL	1725	80024	ea.	24529	

Gas-Tight Syringes for Rheodyne® & Valco® Valves

- Teflon®-tipped plungers.
- Removable needles: point style 3 (90° tip); 2"/51mm.
- Replaceable syringe barrels, plungers, and plunger tips.

Hamilton Syringes

10µL syringe: 26s needle
All others: 22s needle



Volume	Hamilton		
	Model	cat.#	Restek
10µL	1701	21260	
25µL	1702	21261	
50µL	1705	21262	
100µL	1710	21263	
250µL	1725	21264	

SGE Syringes
Needle gauge: 22



Volume	SGE			Restek	
	Model	cat.#	qty.	cat.#	
10µL	10R-GT-LC	24866			
25µL	25R-GT-LC	24867			
50µL	50R-GT-LC	24868			
100µL	100R-GT-LC	24869			
250µL	250R-GT-LC	24870			
500µL	500R-GT-LC	24871			

For many more HPLC and GC syringes, request our catalog or visit our website.