

Replacement Parts for Dionex ASE® Systems

by Neil Mosesman, SPE Product Marketing Manager, and Brad Rightnour, Instrument Innovations Manager

- Designed to meet or exceed performance of original manufacturer's parts.
- Siltek® treated parts assure maximum inertness.
- Save time—order parts when you order GC or HPLC columns and consumables.
- Renowned Restek Plus 1™ service.

Accelerated solvent extraction has become a common technique for fast and reliable extraction of organic materials from solid matrices. Restek now offers several direct replacement parts for ASE® extraction systems, including extraction cells, caps, frits, and PEEK® washers. In addition to stainless steel extraction cells—equivalent to original parts—we offer cells, caps, and frits finished with our innovative Siltek® treatment. Siltek® treatment reduces the surface activity of stainless steel, greatly improving inertness and, therefore, the reliability of analytical results for active compounds. All of our replacement parts for ASE® systems are designed to meet or exceed original equipment manufacturer's performance. Look to Restek to keep your extraction systems running smoothly.

Kits and Parts for ASE® Extraction Units

Each kit includes: 1 extraction chamber, 2 end caps, 2 frits, and 2 PEEK® washers.

Description	Similar to Dionex		
	part #	qty.	cat.#
Extraction Cell Kit for ASE® 200 Extraction Unit, 22mL Tubes	048764	kit	26094
Extraction Cell Kit for ASE® 200 Extraction Unit, Siltek®-Treated, 22mL Tubes	—	kit	26095
Replacement Caps for ASE® 200 Extraction Unit, Universal	049450	2-pk.	26096
Replacement Caps for ASE® 200 Extraction Unit, Siltek®-treated, Universal	—	2-pk.	26097
Extraction Tube for ASE® 200 Extraction Unit, 22mL	048821	ea.	26098
Extraction Tube for ASE® 200 Extraction Unit, Siltek®-treated, 22mL	—	ea.	26099
Replacement Frits for ASE® 200 Extraction Unit, Universal	049453	6-pk.	26100
Replacement Frits for ASE® 200 Extraction Unit, Siltek®-treated, Universal	—	6-pk.	26101


PEEK® Washers for ASE® Extraction Unit

Meet original equipment manufacturer's performance.

Description	Similar to Dionex		
	part #	qty.	cat.#
PEEK® Washers for ASE® 200 Extraction Unit	049454	12-pk.	25256
PEEK® Washers for ASE® 200 Extraction Unit	049454	48-pk.	25257
PEEK® Washers for ASE® 300 Extraction Unit	061687	12-pk.	25393
PEEK® Washers for ASE® 300 Extraction Unit	061687	48-pk.	25394

new!

Siltek® Treated
Parts for ASE®
Systems



Untreated Parts for
ASE® Systems



PEEK® Washers for ASE®
Extraction Unit