

It's Here! The 2005 Restek Catalog!

- 775+ pages / thousands of innovative products.
- Many new chromatograms.
- Helpful technical information.

Some of the new items in the 2005 Restek Catalog:

GC Columns

Rtx[®]-440 - Low-bleed, high-resolution, intermediate-polarity column for many applications.

Rtx[®]-5SII MS, 0.18mm ID - Monitor nanogram levels of semivolatile pollutants in 15 minutes.

Rtx[®]-Dioxin2 - Improved separation of dioxin or furan congeners, compared to diphenyl or high-cyano columns.

Rtx[®]-1PONA Column/Rtx[®]-5PONA Tuning Column - Performance enhanced, for 30% faster SIMDIS analysis, using helium.

Rtx[®]-XLB Column - Ideal for active, higher weight environmental analytes and other compounds.

Too new for our catalog! Read about the Rtx[®]-PCB column on page 13 of this *Advantage*.

HPLC Columns

Allure[™] Aqueous C18 - Excellent retention and selectivity for polar analytes in highly aqueous mobile phases.

Ultra Quat - Monitor paraquat and diquat without ion pairing reagents.

HPLC Method Development Column Kits - Multiple stationary phases help to quickly optimize selectivity.

GC Accessories

Parker Balston FID-1000 Gas Station - Ultra-high purity hydrogen and zero grade air from a single unit.

EZ No-Vent[™] Column/MS Connector for Varian Saturn 2000 Series MSs - Change column in minutes without venting.

Alumaseal[™] aluminum ferrules - The sealing ease of graphite, the security and reliability of metal ferrules.

Vespel[®] ferrules - We give you more choices: Vespel[®], Vespel[®]/graphite, or graphite.

HPLC Accessories

LC Autosampler Syringes

Instrument parts - for Agilent, Beckman, Hitachi, PerkinElmer, Shimadzu, Thermo Separation Products, and Waters equipment.

Mobile Phase Pre-heater - Accurately and reproducibly warms mobile phase, for more consistent analyses.

Restek Performance Coatings

Siltek[®], Sulfiner[®], and Silcosteel[®]-CR treated Swagelok[®] fittings - Inert fittings for demanding applications.

Siltek[®] and Silcosteel[®]-CR treated electropolished stainless steel tubing - For the most inert sample pathways available.

Air Monitoring

Improved SilcoCan[™] and TO-Can[™] air sampling canisters - Superior deactivation; superior protection for the valve.

Canister Heating Jacket - Uniformly heats entire canister, and valve, for faster more efficient cleaning.

Analytical Reference Materials

Calibration mixes, check mixes, surrogates, and internal standards for environmental analyses - volatile organics...

semivolatile organics...pesticides & herbicides...UST monitoring...Canadian environmental methods

Dimethyldichlorosilane (DMDCS) deactivation reagent

Drinking water odor compounds

If you haven't received your copy, request one from your Restek representative.

Or, visit our website: www.restek.com or our new performance coatings website:

www.restekcoatings.com

Restek Trademarks/Service Marks: Allure, Alumaseal, Crossbond, EZ No-Vent, MegaMix, Plus 1, Precision, Press-Tight, Rtx, SeCure, SilcoCan, Silcosteel, Siltek, Sulfiner, TO-Can, Trident, Uniliner, Viva, Vu2 Union, Turning Visions Into Reality, Restek logo. **Other Trademarks:** Aroclor (Monsanto Chemical Co.), ASE (Dionex), Auto SYS (PerkinElmer), Blazer (BASF), Dacthal (Amvac Chemical Corp.), DB (J&W Scientific), PEEK (Victrex plc), Sonar (Sepro Corporation), Swagelok (Crawford Fitting Co.), Terrazole (Uniroyal Chemical Company, Inc.), TRACE (Thermo Finnigan), TrueTube (O'Brien Corporation), Waters (Waters Associates, Inc.), Vespel, Viton, Zorbax (E.I. du Pont de Nemours & Co. Inc.), Vydac (Separations Group, Inc.).

Please direct your comments on this publication to Carrie Sprout at carrie.sprout@restekcorp.com or call 814-353-1300, ext. 2151.



Lit. Cat.# 59077-INT

© 2005 Restek Corporation.

Restek U.S. • 110 Benner Circle • Bellefonte, PA 16823 • 814-353-1300 • 800-356-1688 • fax: 814-353-1309 • www.restek.com

Restek France • phone: 33 (0)1 60 78 32 10 • fax: 33 (0)1 60 78 70 90 • e-mail: restekfr@club-internet.fr

Restek Ireland • phone: 44 2890 814576 • fax: 44 2890 814576 • e-mail: restekeurope@aol.com

Thames Restek U.K. LTD • phone: 44 1494 563377 • fax: 44 1494 564990 • e-mail: sales@thamesrestek.co.uk

Restek GmbH • phone: +49 (0) 6172 2797 0 • fax: +49 (0) 6172 2797 77 • e-mail: info@restekgmbh.de

