

Survival Kits for HPLC

Invaluable for Keeping Your System Running Smoothly!

By Becky Wittrig, Ph.D., HPLC Product Marketing Manager

- Tubing, fittings, and tools for system start-up or maintenance.
- Choose PEEK® or stainless steel components.
- More convenient and more economical than ordering components separately.

Restek HPLC survival kits contain practical selections of tubing, fittings, and tools for setting up or maintaining your HPLC system. The PEEK® Survival Kit contains PEEK® tubing, connectors, and elbows, Teflon® tubing, a tubing cutter and extra blades, a ValvTool wrench, open-end wrenches, and more. The Stainless Steel Survival Kit contains a selection of lengths and IDs of 1/16-inch tubing, plus nuts, ferrules, a ValvTool wrench, and a zero-dead-volume union.

did you know?

Restek offers a wide range of HPLC columns, tools, and accessories, and many replacement parts for Agilent, Beckman, Hitachi, PerkinElmer, Shimadzu, Thermo Separation, and Waters instruments. Call us for a copy of our latest HPLC catalog (lit. cat.# 59241B), or visit us on line.

PEEK® Survival Kit for HPLC

The PEEK® Survival Kit is an invaluable parts kit that contains tools and supplies essential for setting up and maintaining your HPLC system.



Description	qty.	cat.#
PEEK® Survival Kit for HPLC	kit	25322

Stainless Steel Survival Kit for HPLC

Contains a wide range of stainless steel tubing, plus fittings and a ValvTool wrench.



Description	qty.	cat.#
Stainless Steel Survival Kit for HPLC	kit	25097