

Instrument Innovations!

Restek is your #1 source for GC consumables and supplies

by Donna Lidgett, GC Accessories Product Marketing Manager

From Injector to Detector

Restek designs, develops, and markets OEM-equivalent parts and supplies for Agilent, PerkinElmer, Shimadzu, Thermo Finnigan or Varian GC systems. Restek consumables and parts meet or exceed OEM performance, helping you maintain optimum system performance and giving you the convenience and economy of one-stop shopping for all your GC needs.

Check out some of our newest consumables. . .



new! EPC Test Kit for Agilent 6890 GCs

- Kit includes 3 o-rings, 2 plugs, 1 mounting screw and 1 test block.

Description	Similar to Agilent part #	qty.	cat.#
EPC Test Kit for Agilent 6890 GCs	G1530-60960*	kit	24323

*Similar to Agilent part # G1530-60960, but not exact equivalent. Kits differ in parts.



Clear anodized aluminum and high-quality stainless steel.

new! Septum Nut for Shimadzu 17A & 2010 GCs

- One piece design for ease of installation and removal.

Description	Similar to Shimadzu part #	qty.	cat.#
Septum Nut for Shimadzu 17A & 2010 GCs	221-41286-00	ea.	22079
	221-44584-00		



High quality stainless steel.

new! FID Jet for PerkinElmer Auto SYS™ XL

Description	Similar to PE part #	qty.	cat.#
FID Jet for PerkinElmer Auto SYS™ XL	N6100361	ea.	23038



High quality stainless steel.

new! FID Capillary Column Adaptor for PerkinElmer Auto SYS™ XL

Description	Similar to PE part #	qty.	cat.#
For use with PE style capillary nuts			
FID Capillary Column Adaptor for PerkinElmer Auto SYS™ XL	N6120020	ea.	22608
For use with 1/8" compression style nuts			
FID Capillary Column Adaptor for PerkinElmer Auto SYS™ XL	—	ea.	22609



Clear anodized aluminum/
High quality stainless steel.

new! Septum Cap for PerkinElmer Auto SYS™ XL

Description	Similar to PE part #	qty.	cat.#
Septum Cap for PerkinElmer Auto SYS™ XL	N6100153	ea.	22322



High quality stainless steel.

new! Injector Adaptor for PerkinElmer Auto SYS™ XL

Description	Similar to PE part #	qty.	cat.#
For use with PE style capillary nuts			
Injector Adaptor for PerkinElmer Auto SYS™ XL	N6100157	ea.	22318
Siltek®-Treated Injector Adaptor for PerkinElmer Auto SYS™ XL	—	ea.	22320
For use with 1/8" compression style nuts			
Injector Adaptor for PerkinElmer Auto SYS™ XL	—	ea.	22319
Siltek®-Treated Injector Adaptor for PerkinElmer Auto SYS™ XL	—	ea.	22321



Siltek®-treated version for increased inertness.