

# Genuine Restek Supplies & Accessories for ASE® Systems

## Extraction Cell Parts, Collection Vials, PEEK® Washers, Filters

By Neil Mosesman, Sample Preparation Product Marketing Manager

### Extraction Cell Parts for ASE® Systems

- Designed to meet or exceed performance of original manufacturer's parts.
- Polished inner surfaces for easier cleaning; Siltek® deactivation available.

In addition to stainless steel extraction cell parts, we offer bodies, caps, and frits finished with our innovative Siltek® treatment, to greatly improve inertness and, therefore, the reliability of analytical results for active compounds.



Description	Similar to Dionex part #	qty.	Stainless Steel		Siltek®-Treated	
			cat.#	qty.	cat.#	
<b>Parts for ASE® 200 Extraction Cells</b>						
Extraction Cell Body for ASE® 200, 1mL	054973	ea.	26110	ea.	26111	
Extraction Cell Body for ASE® 200, 5mL	054974	ea.	26112	ea.	26113	
Extraction Cell Body for ASE® 200, 11mL	048820	ea.	26114	ea.	26115	
Extraction Cell Body for ASE® 200, 22mL	048821	ea.	26098	ea.	26099	
Extraction Cell Body for ASE® 200, 33mL	048822	ea.	26116	ea.	26117	
Replacement Extraction Cell End Caps for ASE® 200	049450	2-pk.	26096	2-pk.	26097	
Replacement Frits for ASE® 200	049453	10-pk.	26100	10-pk.	26101	
<b>Parts for ASE® 300 Extraction Cells</b>						
Extraction Cell Body for ASE® 300, 10mL		ea.	26172	ea.	26173	
Extraction Cell Body for ASE® 300, 34mL		ea.	26176	ea.	26177	
Extraction Cell Body for ASE® 300, 66mL	056696	ea.	26178	ea.	26179	
Extraction Cell Body for ASE® 300, 100mL	056693	ea.	26132	ea.	26133	
Replacement Extraction Cell End Caps for ASE® 300	056921	2-pk.	26170	2-pk.	26171	
Replacement Frits for ASE® 300/100		6-pk.	26174	6-pk.	26175	

### Accessories for ASE® Systems

Meet original equipment manufacturer's performance.

Description	Similar to Dionex part #	qty.	cat.#
<b>Accessories for ASE® 200 Systems</b>			
PEEK® Washers for ASE® 200	049454	12-pk.	25256
PEEK® Washers for ASE® 200	049454	48-pk.	25257
PEEK® Washers for ASE® 200	049454	250-pk.	26120
Snap Rings for Caps for ASE® 200	049456	10-pk.	26184
Funnel for ASE® 200	056958	ea.	26180
<b>Accessories for ASE® 300/100 Systems</b>			
PEEK® Washers for ASE® 300	061687	12-pk.	25393
PEEK® Washers for ASE® 300	061687	48-pk.	25394
Snap Rings for Caps for ASE® 300/100	056778	12-pk.	26134



PEEK® Washers for ASE® Extraction Units

### 20mm Filters for ASE® 200 Extraction Cells

Consistent porosity, to deliver rapid flow rates and protect the metal frit in the cell from contamination.

Description	Similar to Dionex part #	qty.	cat.#
Cellulose Filters for ASE® 200	049458	100-pk.	26118
Glass Fiber Filters for ASE® 200	047017	100-pk.	26119



### Diatomaceous Earth

Mix with densely packed samples such as clays to improve extraction efficiencies and absorb excess moisture.

Description	Similar to Dionex part #	qty.	cat.#
Diatomaceous Earth, 30/40 mesh	062819	1kg	26033



### 60mL Sample Collection Vials

Cleaned to EPA specifications and supplied assembled with caps and septa.

Description	Similar to Dionex part #	qty.	cat.#
60mL Collection Vials, Clear Glass, for ASE® Systems	048784	72-pk.	26121
60mL Collection Vials, Amber Glass, for ASE® Systems	048781	72-pk.	26122

