

Bulk Adsorbents

For Thorough Sample Preparation and Reliable Results

By Neil Mosesman, Sample Preparation Products Marketing Manager

Florisil® PR

- Pesticide residue grade.
- Each lot certified to meet the requirements of AOAC methodology.
- Packaged in glass containers.

Florisil® PR is commonly used to remove polar interferences from pesticide residues. This bulk material is ideal for labs packing their own chromatography columns for pesticide residue extractions.

Description	qty.	cat.#
Florisil® PR, 60/100 mesh	500gms	26135

Granulated Activated Copper

- Convenient form for removing sulfur from environmental extracts.
- Acidified and activated—ready for use.

Activated copper effectively removes elemental sulfur from environmental extracts. Our acid washed and activated material can be used right out of the package. The 30 mesh granular material eliminates the potential for fine copper particles in filtered extracts.

Description	qty.	cat.#
Granulated Activated Copper, 30 mesh	1kg	26136

Ottawa Sand

- Sample medium for matrix spikes and laboratory control blanks.
- Packaged in convenient 5kg buckets.

Ottawa sand is organics free and is listed in several US EPA methods as the specified medium for matrix spike and laboratory control blanks.

Description	qty.	cat.#
Ottawa Sand	5kg	26137

