

# COOL TOOLS!

Restek Innovations Save You Time and Money

## Spanner Wrench for Agilent 5890/6890 FID Collector Assembly

- Easily remove the nut from the FID collector without damaging the nut.
- Unique, ergonomic handle—easy to grip.



Remove the collector housing...



...easily loosen the nut by aligning the two pins on the bottom of the wrench with the two open slots on the nut...



...then turn counterclockwise...



...and remove.



Description	Similar to Agilent part #	qty.	cat.#
Spanner Wrench for Agilent 5890/6890 FID Collector Assembly	19231-00130	ea.	22329

## Injector Wrench for Agilent 5890/6890/6850 GCs

- Use to remove the septum nut and weldments during GC maintenance.
- High-quality stainless steel construction.
- Meets original equipment performance.



Use the smaller end to remove the septum nut.



Use the larger end to tighten the split/splitless weldment nut.



Description	Similar to Agilent part #	qty.	cat.#
Injector Wrench for Agilent 5890/6890/6850 GCs	19251-00100	ea.	22065

## Injector Wrench for Shimadzu 17A and 2010 GCs

- Designed specifically for removing Shimadzu injection ports.
- High-quality stainless steel construction.



Description	Similar to Shimadzu part #	qty.	cat.#
Injector Wrench for Shimadzu 17A and 2010 GCs	221-46977-00	ea.	21159

## 1/4- to 5/16-inch Open-End Wrench Set

We examined many different wrenches before we decided to offer this high-quality pair for tightening capillary fittings.

Description	qty.	cat.#
Open-End Wrenches (1/4" X 5/16")	2-pk.	20110

