

## Hub-Cap Mobile Phase Accessories

### Simplify Mobile Phase Delivery with the Hub-Cap Filter Kit

Introducing our new Hub Cap filter kit! The Hub-Cap filter allows you to simultaneously transfer and filter your mobile phases. The bottle tops and adaptors are designed to fit securely on 4-liter solvent bottles and eliminate messy, loose-fitting parafilm or foil wraps. Tidy up your mobile phase delivery—try a Hub-Cap today!



**Keep your mobile phase lines under control—use Hub-Cap bottle tops instead of parafilm, aluminum foil, or tape on your mobile phase reservoirs.**

#### Hub-Cap Filter Kit

**new!**

Transfer and filter mobile phases in one easy step.

Kit includes: bottle adapter, bottle adapter nut, filter inlet cap, grid support, vacuum hose barb, tube compression fitting, 47mm grid, 47mm .22 $\mu$ m filter membrane, 47mm .45 $\mu$ m filter membrane, 1/4" OD x 1/8" ID ultra chemical resistant, Teflon® FEP lined Tygon® tubing (3'), 6" x 6" box with shrink wrap insert



cat. #26395

Description	qty.	cat.#
Hub-Cap Filter Kit	kit	26395
Replacement Membrane Filters	qty.	cat.#
Polypropylene Membrane Filters, 47mm, 0.45 $\mu$ m	100-pk.	26396
Polypropylene Membrane Filters, 47mm, 0.22 $\mu$ m	100-pk.	26397
Nylon Membrane Filters, 47mm, 0.45 $\mu$ m	100-pk.	26398
Nylon Membrane Filters, 47mm, 0.22 $\mu$ m	100-pk.	26399

#### Hub-Cap 4 Liter Bottle Tops

Most bottles use a GL45 cap. New Hub-Cap bottle tops are a great way to neatly keep your mobile phase lines where they belong. Use instead of parafilm, aluminum foil, or tape on your mobile phase reservoirs.



cat. #26541

Description	qty.	cat.#
Hub-Cap (assembly of the bottle cap and plug)	kit	26541
Hub-Cap Multi-pack	3-pk.	26542

#### Hub-Cap Adapters

Allow the use of the Opti-Cap™ with 4-liter solvent bottles.



cat. #26538

cat. #26540

Description	qty.	cat.#
Hub-Cap Adapter	ea.	26538
Hub-Cap Adapter Multi-pack	3-pk.	26539
Hub-Cap Adapter and Opti-Cap™	kit	26540