

## NEW! Waste Overflow Indicator for HPLC Systems

By Becky Wittrig, Ph.D, HPLC Product Marketing Manager

- Avoid messy pooling around mobile phase waste containers.
- Audible alarm instantly alerts user, preventing overflow.
- Compact, battery operated unit.

The new Restek Waste Overflow Indicator will help to keep your mobile phase waste where it belongs—in the waste container! Compact, battery operated unit accommodates two lines and fits securely on 4-liter solvent bottles. An audible alarm is given as the solvent waste container approaches capacity, giving you time to empty or change the container. Another innovative design from Restek!

### Waste Overflow Indicator for HPLC Systems

Description	qty.	cat.#
Waste Overflow Indicator for HPLC Systems	kit	26543
Replacement AA Battery for the Waste Overflow Indicator	ea.	26544
Replacement AA Batteries for the Waste Overflow Indicator	3-pk.	26545

For more information about HPLC Accessories, visit us online at [www.restek.com/HPLC](http://www.restek.com/HPLC)

