




Analog Signal Cable

Agilent 5890 GC to Agilent 3396 Integrator
(cat.#s 20650, 20651—replaces Agilent #35900-60610)



Connector 1
front view



Connector 2
front view

Name	Conn. 1 Pin	Conn. 2 Pin	Color	Function
Ground	1	3	Bare	Grounded shield
Negative	2	1	Red	Integrator signal (+)
Positive	3	2	Black	Integrator signal (-)

Note: Connector 1 is not polarized—insert into 5890 socket in either direction.

**Call Technical Service at 800-356-1688 or 814-353-1300, ext. 4
(or your Restek representative) if you have any questions about this
product or any other Restek product.**



#203-03 [007]
Rev. date: 02/08

Restek U.S. • 110 Benner Circle • Bellefonte, PA 16823

phone: 814-353-1300 or 800-356-1688 • fax: 814-353-1309 • www.restek.com

Restek France • phone: 33 (0)1 60 78 32 10 • e-mail: restek@restekfrance.fr

Restek Ireland • phone: 44 2890 814576 • e-mail: resteurope@aol.com

Thames Restek U.K. LTD • phone: 44 1494 563377 • e-mail: sales@thamesrestek.co.uk

Restek GmbH • phone: +49 (0) 6172 2797 0 • e-mail: info@restekgmbh.de