



ALS Stop/Start Cable

(cat.# 20929—replaces Agilent #G1530-60930)



**Connector 1
front view**



**Connector 2
front view**

Conn. 1 Pin	Conn. 2 Pin	Color
1	1	Brown
2	2	Green
3	3	Red
4	4	White
5	5	Black
6	6	Purple
7	7	Yellow
8	8	Orange
9	9	Blue

Note: All pin numbers correspond to pin numbers of other connector.

**Call Technical Service at 800-356-1688 or 814-353-1300, ext. 4
(or your Restek representative) if you have any questions about this
product or any other Restek product.**



#203-03 [025]
Rev. date: 01/07

Restek U.S. • 110 Benner Circle • Bellefonte, PA 16823

phone: 814-353-1300 or 800-356-1688 • fax: 814-353-1309 • www.restek.com

Restek France • phone: 33 (0)1 60 78 32 10 • e-mail: restek@restekfrance.fr

Restek Ireland • phone: 44 2890 814576 • e-mail: resteurope@aol.com

Thames Restek U.K. LTD • phone: 44 1494 563377 • e-mail: sales@thamesrestek.co.uk

Restek GmbH • phone: +49 (0) 6172 2797 0 • e-mail: info@restekgmbh.de