

Personal Sample Inlet

cat.# 22123

Includes:

- 3' x 1/16" OD Teflon® tubing
- Adjustable alligator clip
- Teflon® reducing ferrule
- 1/4" stainless steel nut



Teflon® reducing ferrule allows direct connection of 1/16" tubing and 1/4" flow controller.



Convenient clip can be moved along entire tubing length, for proper attachment in breathing zone.



open



closed



Call Technical Service at 800-356-1688 or 814-353-1300, ext. 4 (or your Restek representative) if you have any questions about this product or any other Restek product.

RESTEK

Chromatography Products
© 2009 Restek Corporation

ISO 9001:2000
Rev. 1.0 FMS0307

#500-13 [002]
Rev. date: 7/09

Restek U.S. • 110 Benner Circle • Bellefonte, PA 16823
phone: 814-353-1300 • 800-356-1688 • fax: 814-353-1309 • www.restek.com

Restek France • phone: +33 (0)1 60 78 32 10 • e-mail: restek@restekfrance.fr

Restek GmbH • phone: +49 (0)6172 2797 0 • e-mail: info@restekgmbh.de

Restek Ireland • phone: +44 (0)2890 814576 • e-mail: restekurope@aol.com

Restek Japan • phone: +81 (3)6459 0025 • e-mail: ryosei.kanaguchi@restek.com

Thames Restek U.K. LTD • phone: +44 (0)1494 563377 • e-mail: sales@thamesrestek.co.uk