

# Hub-Cap Bottle Tops and Adapters

## Mobile Phase Accessories

- Specifically designed for 4-liter bottles!
- Easiest, cleanest way to helium sparge and deliver mobile phases.
- Use Hub-Cap bottle tops instead of parafilm, aluminum foil or tape.



**Chromatography Products**

www.restek.com 800-356-1688 • 814-353-1300

## Hub-Cap Bottle Tops and Adapters


### Hub-Cap 4 Liter Bottle Tops

Description	qty.	cat.#
Hub-Cap (assembly of the bottle cap and plug)	kit	26541
Hub-Cap Multi-pack	3-pk.	26542

### Hub-Cap Adapters

Allows the use of the Opti-Cap™ with 4-liter solvent bottles.

Description	qty.	cat.#
Hub-Cap Adapter	ea.	26538
Hub-Cap Adapter Multi-pack	3-pk.	26539
Hub-Cap Adapter and Opti-Cap™	kit	26540



Use Hub-Cap bottle tops instead of parafilm, aluminum foil or tape on your mobile phase reservoirs.

### Opti-Cap™ Bottle Top

The most economical way to helium-sparge and deliver HPLC mobile phases. The Opti-Cap™ top fits all standard GL-45 bottles and has two 1/8-inch holes and one 1/16-inch hole for tubing.

Description	qty.	cat.#
Opti-Cap™ (Cap and PEEK® Plug)	ea.	25300
Opti-Cap™ Kit (Opti-Cap™, 3 meters of tubing, sparging filters)	kit	25301
Opti-Cap™ Kit with 1L Bottle	kit	25302
Opti-Cap™ Kit with 2L Bottle	kit	25303
<b>Related items and replacement parts</b>	<b>qty.</b>	<b>cat.#</b>
Mobile Phase Sparge Filter: 2µm, stainless steel	ea.	25311
Mobile Phase Inlet Filter: 10µm	ea.	25312
Teflon® Tubing, 1/8" OD x 0.094" ID x 3m (2.4mm ID)	3m	25307
Teflon® Tubing, 1/8" OD x 0.063" ID x 3m (1.6mm ID)	3m	25306
PEEK® Plug, 1/4"-28 threads	3-pk.	25319
1L Graduated Safety-Coated Bottle – GL-45 threads	ea.	25304
2L Graduated Safety-Coated Bottle – GL-45 threads	ea.	25305



Opti-Cap™ Kit with bottle

Restek Trademarks:  
Restek logo.

Other Trademarks:  
Opti-Cap (Jour Research), PEEK (Victrex plc), Teflon (E.I. du Pont de Nemours & Co., Inc.)



**Restek U.S.** • 110 Benner Circle • Bellefonte, PA 16823 • 814-353-1300 • 800-356-1688 • fax: 814-353-1309 • [www.restek.com](http://www.restek.com)  
**Restek France** • phone: 33 (0)1 60 78 32 10 • fax: 33 (0)1 60 78 70 90 • e-mail: [restek@restekfrance.fr](mailto:restek@restekfrance.fr)  
**Restek Ireland** • phone: 44 2890 814576 • fax: 44 2890 814576 • e-mail: [restekeurope@aol.com](mailto:restekeurope@aol.com)  
**Thames Restek U.K. LTD** • phone: 44 1494 563377 • fax: 44 1494 564990 • e-mail: [sales@thamesrestek.co.uk](mailto:sales@thamesrestek.co.uk)  
**Restek GmbH** • phone: +49 (0) 6172 2797 0 • fax: +49 (0) 6172 2797 77 • e-mail: [info@restekgmbh.de](mailto:info@restekgmbh.de)

Lit. Cat.# 580144

© 2007 Restek Corporation.

